



Dragon Crown Medical Co., Ltd.


Tel.: +0086-531-51217461

Web.: www.dcvertebrae.com

Add.: No.978, Tianchen Road,
Hi-Tech Development Zone,
Jinan, Shandong Province, China.




RF Electrodes



Dragon Crown Medical Co., Ltd.

General Description




Model	DJ - I	DJ - II	DJ - III	DJ - IV
Working length	L=350mm	L=320mm	L=230mm	L=200mm
OD	φ2.5mm	φ2.5mm	φ2.5mm	φ2.5mm

Product Features

- RF Electrodes connect to a RF Generator. Used to cut living tissue or perform coagulation or ablation.
- The product is compatible with most RF Generators. The handle design is more ergonomic, easy to use, flexible and precise.
- During operation, the relative position between the two poles remains unchanged, making the ablation discharge more stable.
- Adapter plugs can be customized according to different brands of RF Generators.
- The product has been upgraded many times. The head is heat-processed and has high elasticity. The electrode pads are ground with high precision . Not causing soft tissue adhesion.
- The product is sterily packaged and is for one-time use.




RF Generator



Product Parameters

- Parameters of RF Generator used together: Rated accessory voltage: 750V, corresponding gear output power: ≤ 120W
- The dimensions of the working part of RF Electrodes are shown in the figure.
- Working diameter 2.5mm



Suitable for DragonCrown surgical instruments



With PELD Surgical System



With Spinal Endoscopic Large Channel Surgical System



With Cervical Spine Surgical System



With DMSE System



With VBE System



With UBE System